



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY Speaking

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

FITC Adsorbed Rabbit anti-Mouse Macrophage

CLFAD31240

Lot #: 6304F

Description: Cedarlane's **CLFAD31240** is a FITC-conjugated rabbit IgG fraction directed against mouse macrophages. This fraction was produced from an antiserum that had been adsorbed repeatedly with mouse red blood cells, thymocytes and spleen cells to remove most antibodies to antigens common among different cell types. This adsorption significantly improves specificity. The FITC conjugate has been purified by removing free fluorochrome, under-conjugated and over-conjugated antibodies. It has been analyzed for fluorochrome:protein ratio.

Presentation: 2 mg in sterile phosphate buffered saline.

Protein Concentration: 2.5 mg/ml

Fluorescein:Protein Ratio: 2.1:1

Storage: Divide into useful aliquots and store at -20°C. Avoid prolonged exposure to light.

Preservatives: 0.05% sodium azide

Applications: This antibody is suitable for use in flow cytometry, immunofluorescence and cryostat sections.

Specificity: This FITC antiserum fluoresces mouse macrophages but not mouse thymocytes or splenocytes.

Sensitivity: This lot (6304F) of antiserum fluoresces at dilutions ranging from neat to 1:30.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

LW 05/10/01

For Research Use Only

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com