



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

FITC Adsorbed Rabbit anti-Mouse Thrombocyte

CLFAD31440

Lot #: 6306AF-F

DESCRIPTION:

Cedarlane's **CLFAD31440** is a FITC-conjugated rabbit IgG fraction directed against mouse Thrombocytes. This fraction was produced from an antiserum that had been adsorbed repeatedly with mouse red blood cells, thymocytes and spleen cells to remove most antibodies to antigens common among different cell types. This adsorption significantly improves specificity. The FITC conjugate has been purified by removing free fluorochrome, under-conjugated and over-conjugated antibodies.

PRESENTATION:

2 mg in sterile phosphate buffered saline.

PROTEIN CONCENTRATION:

2.0 mg/ml

FLUORESEIN:PROTEIN RATIO:

2.1 mole to 1.0 mole

STORAGE:

Divide into useful aliquots and store at -20°C. Avoid prolonged exposure to light. For frequent use store at +4°C.

PRESERVATIVES:

0.05% sodium azide

APPLICATIONS:

This antibody is suitable for use in flow cytometry, immunofluorescence and cryostat sections.

SPECIFICITY:

The above lot of conjugate fluoresces mouse thrombocytes but does not fluoresce chicken RBCs.

SENSITIVITY:

The above lot of conjugate fluoresces at dilutions 1:1 to 1:250

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Laboratory Research For Research Use Only

MW 09/30/15

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com