



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.  
Please contact CEDARLANE® for lot specific information.

## FITC Rabbit anti-Swine Albumin

**CLFAG16140**

**Lot:** XXXXX

**Description:** Cedarlane's CLFAG16140 is prepared by the covalent coupling of fluorescein to the IgG fraction of high titered rabbit antiserum. The conjugates are purified to remove free fluorochrome, under-conjugated and over-conjugated antibodies. The fluorescent products are analyzed for their fluorochrome: protein ratios.

**Fluorescein:Protein Ratio:** TBD

**Presentation:** 3 mg of antibody in phosphate buffered saline containing 0.05% sodium azide as preservative.

**Protein Concentration:** TBD

**Storage:** Avoid prolonged exposure to light. Stable at 4°C. For long term storage, divide into useful aliquots and store at -20°C.

**Applications:** This FITC conjugated antibody can be used in flow cytometry and immunofluorescence.

### Ouchterlony: Results

|                         | Test                            |
|-------------------------|---------------------------------|
| Sample Size:            | 10 µl                           |
| Distance Between Wells: | 5 mm                            |
| Test time:              | 24 hours                        |
| Sample tested:          | swine albumin, purified protein |
| Results:                | TBD                             |

**Laboratory Reagent For Research Use Only**

JV 05/30/07

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
*in North America*

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)