



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

FITC Rabbit anti-Human Albumin

CLFAG2140
Lot: 2001G-2F

DESCRIPTION:

Cedarlane's CLFAG2140 is prepared by the covalent coupling of fluorescein to the IgG fraction of high titered rabbit antiserum. The conjugates are purified to remove free fluorochrome, under-conjugated and over-conjugated antibodies. The fluorescent products are analyzed for their fluorochrome: protein ratios.

FLUORESCHEIN:PROTEIN RATIO:

2.1:1.0 mol

PRESENTATION:

3 mg of antibody in phosphate buffered saline containing 0.05% sodium azide as preservative.

PROTEIN CONCENTRATION:

3.1 mg/ml.

STORAGE:

Avoid prolonged exposure to light. Stable at 4°C. For long term storage, divide into useful aliquots and store at -20°C.

APPLICATIONS:

This FITC conjugated antibody can be used in flow cytometry and immunofluorescence.

OUCHTERLONY: RESULTS

	Test
Sample Size:	10 µl
Distance Between Wells:	5 mm
Test time:	24 hours
Sample tested:	human albumin, purified protein
Results:	one fluorescent band

For Research Use Only

JV 04/05/10

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com