



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

FITC Rabbit anti-Human IgG

CLFAG2340
Lot #: 0135SG-HF

TEST TYPE:

Fluorescein: Protein ratio = 2.3mol FITC: 1mol IgG.

Protein Concentration: 0.3 mg IgG/ml in PBS.

Agglutination: Agglutinates and fluoresces human IgG tanned fixed sheep RBC's.

PRESENTATION:

3 mg of frozen antibody in phosphate buffered saline containing 0.05% sodium azide as preservative.

STERILITY TEST:

Thioglycollate Broth @ 37°C
Trypticase Soy Broth @ 37°C
Sabourad's Broth @ 22°C
Results: No Test

STORAGE:

For long term storage: Divide into usable aliquots and store at -20°C. For frequent use: Store at 4-8°C.

For Research Use Only

BA 01/11/12

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com