



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.
Please contact CEDARLANE® for lot specific information.

FITC Rabbit anti-Mouse Transferrin

CLFAG3240

Lot #: 3033F-M

Description:

Cedarlane's rabbit anti-mouse transferrin is prepared by the covalent coupling of fluorescein to the IgG fraction of high titered rabbit antiserum. The conjugates are purified to remove free fluorochrome, under-conjugated and over-conjugated antibodies.

Fluorescien: Protein Ratio: 1.4:1

Presentation:

3 mg of antibody in phosphate buffered saline (pH 7.3) containing 0.05% sodium azide (NaN₃) as preservative.

Protein Concentration: 3.0 mg/ml.

Storage:

Avoid prolonged exposure to light. Stable at 4°C. For long term storage, divide into useful aliquots and store at -20°C.

Applications:

This FITC conjugated antibody can be used in flow cytometry and immunofluorescence.

Ouchterlony: Results

	Test #1
Sample Size:	10 µl
Distance Between Wells:	10 µl
Test time:	24 hours
Sample tested:	mouse transferrin, purified protein
Results:	one fluorescent band

For Research Use Only

JV 04/13/07

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com