



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## FITC Rabbit anti-Rat Albumin

**CLFAG5140**

**Lot:** 3006F

### DESCRIPTION:

Cedarlane's CLFAG5140 is prepared by the covalent coupling of fluorescein to the IgG fraction of high titered rabbit antiserum. The conjugates are purified to remove free fluorochrome, under-conjugated and over-conjugated antibodies. The fluorescent products are analyzed for their fluorochrome: protein ratios.

### FLUORESCHEIN:PROTEIN RATIO:

1.64:1

### PRESENTATION:

3 mg of antibody in phosphate buffered saline (pH 7.3) containing 0.05% sodium azide as preservative.

### PROTEIN CONCENTRATION:

5 mg/ml.

### STORAGE:

Avoid prolonged exposure to light. Stable at 4°C. For long term storage, divide into useful aliquots and store at -20°C.

### APPLICATIONS:

This FITC conjugated antibody can be used in flow cytometry and immunofluorescence.

### OUCHTERLONY: RESULTS

	<b>Test #1</b>	<b>Test #2</b>
Sample Size:	10 µl	10 µl
Distance Between Wells:	5 mm	5 mm
Test time:	24 hours	24 hours
Sample tested:	rat albumin, purified protein	normal rat serum
Results:	one fluorescent band	one band homogeneous with band in test #1

**Laboratory Reagent For Research Use Only**

JV 10/17/07

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)