



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Human Complement Component C1q

CLPRO554L
CLPRO554L-2
CLPRO554L-3

Description: Human C1q produced in Human plasma having a molecular mass of 410 kDa. The C1 component of complement is a complex of 2 C1q sub-units and a tetramer comprising 2 C1r and 2 C1s sub-units. C1 activation is regulated by an inhibitor (C1 inhibitor). Reductions in the levels of components C3, C4, C1 and C2 can be associated with glomerular disease, endocarditis and systemic lupus erythematosus (SLE).

Synonyms: Component C1q, Complement C1q, Complement Component C1q, C1q

Source: Human Plasma

Presentation: 200 µg (**CLPRO554L**), 1 mg (**CLPRO554L-2**), or 10 mg (**CLPRO554L-3**) sterile filtered solution (1 mg/ml). 10 mM HEPES and 300mM NaCl, pH 7.2.

Stability: Human C1q is stable at 4°C if entire vial will be used within 2-4 weeks. Store, frozen below -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Purity: Greater than 95.0% as determined by SDS-PAGE.

Human Virus Test: Plasma from each donor has been tested and found negative for antibody to HIV-1, HIV-2, HCV, HBSAG.

Laboratory Reagent For Research Use Only

MW 01/26/23

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com