



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Human Complement Component C4c

CLPRO556
CLPRO556-2
CLPRO556-3

Introduction: Complement C4c is a degradation product of C4b which is cleaved by C4b/C3b inactivator to yield C4d and C4c. Complement C4c interacts with C1 and C2 to form C3 convertase of the classic activation pathway. Systemic lupus erythematosus (SLE) is often associated with congenital C4 deficiency. Reduced levels of components of the classic pathway (C1, C2, C4, C3) are common after such activation, e.g. in SLE, acute serum sickness and conditions associated with circulating immune complexes.

Description: Human C4c produced in Human plasma having a molecular mass of 210 kDa.

Synonyms: Complement C4c, Complement Component C4c, C4c.

Source: Human Plasma.

Presentation: 200 µg (CLPRO556), 1 mg (CLPRO556-2), or 10 mg (CLPRO556-3), sterile filtered clear solution. The Human Complement C4c solution is in Sodium Phosphate buffer pH 7.0 containing 0.15M NaCl and 0.09% NaN₃.

Stability: Human C4c although stable at 4°C for 1 week, should be stored at -15°C. Please avoid freeze-thaw cycles.

Purity: Greater than 99.0%.

Human Virus Test: Plasma from each donor has been tested and found negative for HIV-1 & 2 antibodies, Hepatitis B surface antigen, and Hepatitis C antibodies.

Laboratory Reagent For Research Use Only

JK 10/27/21

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com