



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Texas Red Rabbit anti-Rat PMN

CLRAD51140

Lot: 5845AD-TR

DESCRIPTION:

Cedarlane's CLRAD51140 is prepared by the covalent coupling of Texas Red to the IgG fraction of high titered rabbit antiserum. The conjugates are purified to remove free fluorochrome, under-conjugated and over-conjugated antibodies. The fluorescent products are analyzed for their fluorochrome: protein ratios.

TEXAS RED:PROTEIN RATIO:

(Absorbance ratio 596:280 nm) = 0.69 moles:1.0 mole

PRESENTATION:

2 mg of antibody in phosphate buffered saline containing 0.05% sodium azide as preservative.

PROTEIN CONCENTRATION:

2.1 mg/ml

STORAGE:

Avoid prolonged exposure to light. Stable at 4°C. For long term storage, divide into useful aliquots and store at -20°C.

APPLICATIONS:

This Texas Red conjugated antibody can be used in flow cytometry and immunofluorescence.

SPECIFICITY:

The above lot of conjugate fluoresces and agglutinates fresh rat PMNs but does not fluoresce or agglutinate fixed chicken RBCs.

SENSITIVITY:

The above lot of conjugate fluoresces at dilutions neat to 1:20

For Research Use Only

MW 11/18/15

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com