



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®]
www.cedarlanelabs.com



Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Texas Red Rabbit anti-Rat Albumin

CLRAG5140

Lot: 3006TR

DESCRIPTION:

Cedarlane's CLRAG5140 is prepared by the covalent coupling of Texas Red to the IgG fraction of high titered rabbit antiserum. The conjugates are purified to remove free fluorochrome, under-conjugated and over-conjugated antibodies. The fluorescent products are analyzed for their fluorochrome: protein ratios.

TEXAS RED:PROTEIN RATIO:

(Absorbance ratio 596:280 nm) = 0.2 moles:1.0 mole

PRESENTATION:

3 mg of antibody in phosphate buffered saline containing 0.05% sodium azide as preservative.

PROTEIN CONCENTRATION:

3 mg/ml

STORAGE:

Avoid prolonged exposure to light. Stable at 4°C. For long term storage, divide into useful aliquots and store at -20°C.

APPLICATIONS:

This Texas Red conjugated antibody can be used in flow cytometry and immunofluorescence.

OUCHTERLONY: RESULTS

	Test 1
Sample Size:	10 µl
Distance Between Wells:	5 mm
Test time:	24 hours
Sample tested:	rat albumin
Results:	positive red fluorescent band

For Research Use Only

MW 03/01/16

Visit our website for your local distributor.

CEDARLANE[®]



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com