



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Normal Rabbit Serum, Mouse Cell Adsorbed

CLSD403-M
Lot: 0924MAD

DESCRIPTION:

Rabbit serum is prepared from normal, non-immunized rabbits. It is processed from whole blood. This normal rabbit serum is adsorbed with mouse red blood cells, thymocytes, and spleen cells using the same procedures done to make our adsorbed rabbit antisera to mouse macrophage, neutrophil, and thrombocyte. This adsorbed normal rabbit serum is the control serum when using the above three adsorbed antisera in immunological tests in vivo or in vitro.

PRESENTATION:

5 ml, clear liquid with a slight straw-tan colour.

STORAGE:

Store at -20°C.

SHELF-LIFE:

When stored properly, this product will keep its usefulness for 2 years.

PRESERVATIVES:

None (0.2 µm filtered)

STERILITY TEST:

None

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Laboratory Reagent For Research Use Only

MW 05/28/15

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com