



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Normal Rabbit Serum, Mouse Cell Adsorbed, IgG Fraction

**CLSDG403-M**

**Lot:** 0920MAD-GG

**Description:** Rabbit serum is prepared from normal, non-immunized rabbits. It is processed from whole blood. This normal rabbit serum is first adsorbed with mouse red blood cells, thymocytes, and spleen cells. The IgG fraction is then made from this adsorbed normal rabbit serum. These are the same procedures used to make our IgG fractions of adsorbed rabbit antisera to mouse macrophage, neutrophil, and thrombocyte cells. The IgG fractions of adsorbed normal rabbit serum is the control serum when using the above three IgG fractions of adsorbed antisera in immunological tests in vivo or in vitro.

**Presentation:** 20 mg (5.0 mg IgG/ml in sterile PBS).

**Storage:** Store at -20°C.

**Shelf-Life:** When stored properly, this product will keep its usefulness for 2 years.

**Preservatives:** None (0.2 µm filtered)

**Sterility test:** None

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

**Laboratory Reagent For Research Use Only**

MW 06/15/16

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)