



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY Speaking

Place your order with CEDARLANE® or your local distributor.  
Please contact CEDARLANE® for lot specific information.

## Mouse Anti-Human CA 125 Monoclonal Antibody

CLSG36531  
CLSG36531-05  
Lot:

<b>Immunogen:</b>	Partially purified human mucin fraction from a pool of cancer tissues of patients with epithelial ovarian.
<b>Clone:</b>	OV185:1
<b>Isotype:</b>	IgG1
<b>Format:</b>	1ml (CLSG36531) or 0.5 ml (CLSG36531-05). This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.
<b>Staining Protocol:</b>	We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. <u>Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining.</u>
<b>Specificity:</b>	This antibody reacts specifically with CA125 ovarian cancer antigen. Several studies have shown that CA125 is useful tumor marker for ovarian epithelial malignancies. However, CA125 has been described in other neoplasms such as seminal vesicle and anaplastic lymphomas.
<b>Positive Control:</b>	Ovarian carcinoma
<b>Cellular Localization:</b>	Extracellular membrane
<b>Storage:</b>	Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

*Continued...*

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
**in North America**

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.  
Please contact CEDARLANE® for lot specific information.

## References:

- i) Alagoz et al. Gynecol Oncol 53:93, 1994.
- ii) Franchi et al. Italian J Gynecol Obstet 5: 149, 1993.
- iii) Marci et al. Gynecol Obstet Invest 37: 143, 1994.
- iv) Ohmori et al. Pathol Intl 4: 333, 1994.

JV 10/06/08

**Laboratory Reagent For Research Use Only**

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



*toll free: 1-800-268-5058*  
*in North America*

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)