



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Mouse Anti-CD34 Endothelial Cell Monoclonal Antibody

CLSG36584
CLSG36584-05
Lot:

Immunogen: Detergent solubilized vesicular suspension prepared from a perfusate of human term placenta.

Clone: QBEND/10

Isotype: IgG1, Kappa

Format: 1 ml (**CLSG36580**) or 0.5 ml (**CLSG36580-05**). This antibody is supplied as purified immunoglobulin containing sodium azide as a preservative.

Titer/Working Dilution: This antibody may be diluted to a titer of 1:25-1:100 in an ABC method. The final dilution should be determined by the user based upon the staining conditions employed.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. Suitable for formalin fixed, paraffin embedded tissue sections.

Specificity: This antibody reacts with a glycoprotein of MW 110 kD, present on the surface of hematopoietic cells (i.e. myeloid and lymphoid progenitors). It stains the vascular endothelial cells.

Positive Control: Tonsil

Cellular Localization: Cell membrane

Continued...

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com

Storage:

Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

References:

- i) Anthony et al. J Clin Pathol 44: 29, 1991.
- ii) Fina et al. Blood 75: 2417, 1990.
- iii) Ramani et al. Histopathol 17: 237, 1990.

Laboratory Reagent For Research Use Only

BA 11/06/12