



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Mouse Anti-CD34 Endothelial Cell Monoclonal Antibody

**CLSG36584**  
**CLSG36584-05**  
**Lot:**

**Immunogen:** Detergent solubilized vesicular suspension prepared from a perfusate of human term placenta.

**Clone:** QBEND/10

**Isotype:** IgG1, Kappa

**Format:** 1 ml (**CLSG36580**) or 0.5 ml (**CLSG36580-05**). This antibody is supplied as purified immunoglobulin containing sodium azide as a preservative.

**Titer/Working Dilution:** This antibody may be diluted to a titer of 1:25-1:100 in an ABC method. The final dilution should be determined by the user based upon the staining conditions employed.

**Staining Protocol:** We suggest an incubation period of 30 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. Suitable for formalin fixed, paraffin embedded tissue sections.

**Specificity:** This antibody reacts with a glycoprotein of MW 110 kD, present on the surface of hematopoietic cells (i.e. myeloid and lymphoid progenitors). It stains the vascular endothelial cells.

**Positive Control:** Tonsil

**Cellular Localization:** Cell membrane

*Continued...*

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)

**Storage:**

Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

**References:**

- i) Anthony et al. J Clin Pathol 44: 29, 1991.
- ii) Fina et al. Blood 75: 2417, 1990.
- iii) Ramani et al. Histopathol 17: 237, 1990.

**Laboratory Reagent For Research Use Only**

BA 11/06/12