



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

*Please contact CEDARLANE® for lot specific information.*

## Purified Mouse Anti-Human Desmin

**CLSG36697**

**Lot:**

**Immunogen:**

BALB/C mice were immunized with purified desmin from human muscle.

**Clone:**

D33

**Isotype:**

IgG1, Kappa

**Format:**

This antibody is supplied as tissue culture supernatant containing sodium azide as preservative.

**Titer/Working Dilution:**

This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.

**Staining Protocol:**

We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed,

*Continued Overleaf...*

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



*toll free: 1-800-268-5058*  
*in North America*

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)

optimal incubation should be determined by the user. Suitable for formalin/paraffin tissues.

- Specificity:** This antibody reacts with desmin protein of 53 kD. Both striated as well as smooth muscle cells stain well with this antibody. The antibody shows a broad interspecific cross-reactivity i.e. reacts with human, rat, chicken and hamster.
- Positive Control:** Leiomyoma
- Cellular Localization:** Cytoplasmic
- Other applications:** Western Blotting (1:50 to 1:100)
- Storage:** Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.
- References:**
- i) Altmannsberger et al. Am J Pathol 118: 85, 1985.
  - ii) Mieltien et al. Arch Dermatol 121: 736, 1985.
  - iii) Gatter et al. J Clin Pathol 39: 950, 1986.

**For Research Use Only. Not for Diagnostic or Therapeutic Work.**