



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY Speaking

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Mouse Anti-Human Glycophorin C

CLSG36728

Lot:

Description:

This monoclonal antibody is specific for human beta sialoglycophorin (Glycophorin C). Glycophorin C is found on the human erythrocyte membrane and is one of four sialic acid rich polypeptides.

Presentation:

1.0 ml, tissue culture supernatant containing 0.05% sodium azide (NaN_3) as a preservative.

Immunogen:

BALB/C mice were injected with Red cell ghosts.

Isotype:

Mouse IgG₁ kappa

Clone:

Ret40f

Applications:

This antibody is suitable for use in immunohistochemistry (formalin/paraffin tissues). We recommend a starting dilution of 1:15 – 1:30 for use with an ABC method, with an incubation period of 30 minutes at room temperature.

For optimal results in various applications, it is recommended that each investigator determine dilutions and incubation period appropriate for individual use.

Positive Control:

Tonsil

Cellular Localization:

Cell Membrane

Storage:

4°C. Check label for expiry date.

References:

1. Gatter *et al.* Histopathol. 13: 257, 1988.

FOR RESEARCH USE ONLY

TN 05/31/05

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com