



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Mouse Anti-Human Laminin Monoclonal Antibody

**CLSG36765**  
**CLSG36765-05**  
**Lot:**

**Immunogen:** The A chain of human laminin.

**Clone:** 4C7

**Isotype:** IgG2a, kappa

**Format:** 1 ml (**CLSG36765**) or 0.5 ml (**CLSG36765-05**). This antibody is supplied as tissue culture supernatant fraction containing sodium azide as preservative.

**Titer/Working Dilution:** This antibody may be diluted to a titer of 1:5-1:20 in an ABC method. The final dilution should be determined by the user based upon the staining conditions employed.

**Staining Protocol:** We suggest an incubation period of 30-60 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. Proteolytic treatment of paraffin tissue sections is recommended.

**Specificity:** This antibody reacts with the terminal globular domain of the A chain. It does not cross react with Laminin from mouse or rat. The staining with this antibody is sensitive to fixation of paraffin tissues and use of Proteolytic enzymatic treatment.

**Positive Control:** Skin

**Cellular Localization:** Basement membrane

*Continued...*

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)

**Storage:**

Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

**References:**

- i) Autio-Harminen et al. Lab Invest 64L 483, 1991.
- ii) Engval et al. Cell Regulation 1: 731, 1990.
- iii) Kleinman et al. Proc Natl Acad Sci USA 85: 1282, 1988.

BA 03/27/12

**Laboratory Reagent For Research Use Only**