



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

## Purified Mouse Anti-Porcine Vimentin Monoclonal Antibody

CLSG36877

CLSG36877-05

Lot:

**Immunogen:**

BALB/C mice were injected with purified vimentin from porcine eye lens.

**Clone:**

V9

**Isotype:**

IgG1, Kappa

**Format:**

1 ml (CLSG36877) or 0.5 ml (CLSG36877-05). This antibody is supplied as purified immunoglobulin from the ascites by Protein G chromatography and contains sodium azide as preservative.

**Titer/Working Dilution:**

This antibody may be diluted to a titer of 1:25-1:50 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.

**Staining Protocol:**

We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining.

**Specificity:**

This antibody reacts with the 57 kDa intermediate filament protein present in the cells of mesenchymal origin. The antibody does not show any cross-reactivity with other intermediate filament proteins. However, it shows a broad interspecies cross-reactivity.

**Positive Control:**

Sarcoma

**Cellular Localization:**

Cytoplasm

**Other applications:**

Immunoprecipitation (1:50), Western Blotting (1:50 to 1:100)

**Storage:**

Store the antibody at 2-8°C. Do not use it beyond the expiration date stated on the label.

**References:**

- i) Osborn et al. J Cell Biol 34: 137, 1984.
- ii) Azumi et al. Am J Pathol 88: 286, 1987
- iii) Buley et al. J Clin Pathol 40: 136, 1987.

**Laboratory Reagent For Research Use Only**

JV 04/03/08

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**

**in North America**

**phone: (905) 878-8891 • fax: (905) 878-7800**

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

**website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)**