



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY Speaking

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Rabbit Anti-*Helicobacter Pylori*

CLSG36942

Lot #:

Immunogen: Total lysate of *H. pylori* were injected to the rabbit.

Host: Rabbit

Size: 1 ml

Format: Purified immunoglobulin fraction of rabbit antiserum against human *H. pylori* containing sodium azide as preservative.

Titer/Working Dilution: This antibody may be diluted to a titer of 1:50-1:100 in an ABC method. However, depending on the staining conditions employed, we suggest that final dilution should be determined by the user.

Staining Protocol: We suggest an incubation period of 30 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation conditions should be determined by the user. Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining.

Specificity: This antibody reacts with the whole *H. pylori*. The antibody stains the bacteria present on the surface of the epithelium or in the cytoplasm of the epithelial cells. *H. pylori* is known to cause peptic ulcer and chronic gastritis in human.

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com

Positive Control: Infected stomach

Cellular Localization: Cytoplasmic

Storage: Store the antibody at 4-8°C. Do not use it beyond the expiration date stated on the label.

References: i) Andersen *et al.* *Apmis* **96**: 559, 1988.

For Research Use Only. Not for Diagnostic or Therapeutic Work.

CH 09/02/05