



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>  
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)



Delivering You Today's Innovations for  
the Science of Tomorrow<sup>™</sup>

## Rabbit Anti-Human Thyroglobulin Polyclonal Antibody

**CLSG36992**  
**CLSG36992-05**  
**Lot:**

**Immunogen:** Thyroglobulin isolated from human thyroid glands.

**Host:** Rabbit

**Format:** 1 ml (**CLSG36992**), or 0.5 ml (**CLSG36992-05**). Purified immunoglobulin fraction of rabbit antiserum against human thyroglobulin containing sodium azide as a preservative.

**Titer/Working Dilution:** This antibody may be diluted to a titer of 1:50-1:200 in an ABC method. The final dilution should be determined by the user based upon the staining conditions employed.

**Staining Protocol:** We suggest an incubation period of 30 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. Suitable for formalin fixed, paraffin embedded tissue sections.

**Specificity:** This antibody reacts with human thyroglobulin containing cells in thyroid.

**Positive Control:** Thyroid

**Cellular Localization:** Cytoplasmic

**Storage:** Store at 2-8°C. Do not use beyond the expiration date stated on the label.

**References:**  
i) Burt et al. J Clin Pathol 38:893, 1985.  
ii) Lauberg et al. Scand J Clin Lab Invest 47: 685, 1987.  
iii)Kastrup et al. Scand J Clin Lab Invest 45: 471, 1985.

**Laboratory Reagent For Research Use Only**

JV 05/10/11

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup>



[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)