



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>  
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)



Delivering You Today's Innovations for  
the Science of Tomorrow™

## Rabbit anti-Progesterone Receptor Polyclonal Antibody

**CLSG80010**  
**CLSG80010-05**  
**Lot:**

- Immunogen:** A synthetic peptide derived from the C-terminal of human progesterone receptor.
- Host:** Rabbit
- Format:** 1 ml (**CLSG80010**), or 0.5 ml (**CLSG80010-05**). Purified immunoglobulin fraction of rabbit antiserum against PgR protein containing sodium azide as a preservative.
- Specificity:** Human PgR exists in two forms: 116kDa (B-form) and 81kDa (A-form). It acts as a ligand activated transcription factor to regulate expression of target genes. Null mutation in the PgR gene leads to pleiotrophic reproductive abnormalities.
- Positive Control:** Breast carcinoma
- Cellular Localization:** Nuclear
- Storage:** Store at 2-8°C. Do not use beyond the expiration date stated on the label.

**Laboratory Reagent For Research Use Only**

JV 09/07/10

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup>



[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)