



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.  
*Please contact CEDARLANE® for lot specific information.*

## Purified Mouse Anti-Human Daxx Monoclonal Antibody

**CLX182AP**  
**LOT:**

**Size:** 0.1 mg

**Clone:** DAXX-01

**Isotype:** Mouse IgG1

**Specificity:** Daxx.

**Immunogen:** Recombinant C-terminal part (aa 558-740) of human Daxx.

**Species Reactivity:** Human

**Application:** Western Blotting, Immunocytochemistry.

**Purity:** > 95% by SDS-PAGE

**Purification:** Purified from ascites by protein-A affinity chromatography.

**Concentration:** 1 mg/ml

**Storage Buffer:** Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4

**Storage / Stability:** Store at 2-8°C. Do not use after expiration date stamped on vial label. For long-term storage, aliquot and store at -20°C. Avoid freeze/thaw cycles.

**Expiration:** See vial label

**Lot Number:** See vial label

**References:** Unpublished.

*This product is made available in collaboration with Exbio.*

**For laboratory research only, not for drug, diagnostic or other use**

JCr 05/25/06

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
**in North America**

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
**website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)**