



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.  
Please contact CEDARLANE® for lot specific information.

## Purified Mouse Anti-Human CD36 Monoclonal Antibody

**CLX238AP**

**Lot:**

**Size:** 0.1mg

**Clone:** TR9

**Isotype:** Mouse IgG1

**Specificity:** The antibody TR9 reacts with CD36 (GPIIb), an 85 kDa integral membrane glycoprotein expressed on platelets, macrophages, endothelial cells, early erythroid cells and megakaryocytes. This antibody cross-blocks binding of FITC-labelled standard antibody OKM5.

CD36 antibodies inhibit adhesive functions (e.g. adherence of infected erythrocytes to target cells).

**Immunogen:** Platelets

**Species Reactivity:** Human

**Application:** Flow Cytometry

**Purity:** > 95% (by SDS-PAGE)

**Purification:** Purified from ascites by precipitation methods and ion exchange chromatography.

**Concentration:** 1 mg/ml

**Storage Buffer:** 20 mM Tris/HCl with 15 mM sodium azide, approx. pH 7.4

**Storage / Stability:** Store at 2-8°C. Do not use after expiration date stamped on vial label. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.

*This product is made available in collaboration with Exbio Praha, a.s.*

**For Research Use Only**

JCr 9/12/06

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
**in North America**

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information  
**website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)**