



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Monoclonal Antibody to Blood Group A Hybridoma Culture Supernatant

CLX255

Size: 1.0 mL

Clone: HE-195

Isotype: Mouse IgM

Specificity: The antibody HE-195 recognizes human blood group A including weak variants A₃, A_X, A₃B, A_XB. The specificity of antibody HE-195 was confirmed by comparison of specificity and reactivity to standard reagent using >5.000 samples of blood.

Immunogen: Mixture of erythrocytes of group A1 and glycoprotein fraction isolated from saliva of secretors with blood group A.

Species Reactivity: Human

Application: **Agglutination**
Working concentrations should be determined by the investigator.

Preparation: The hybridoma culture supernatant is 4x concentrated by ultrafiltration using 100 kDa-cut off membrane.

Storage / Stability: Store at 2-8°C. Do not use after expiration date stamped on vial label. Do not freeze.

Background: Human blood group A antigen belongs to a group of carbohydrate determinants carried on both glycolipids and glycoproteins; it is detected on erythrocytes and certain epithelial cells.

References: *Nemec M, Drimalova D, Horejsi V, Vanak J, Bartek J, Viklicky V: Murine monoclonal antibodies to human A erythrocytes: differential reactivity with N-acetyl-D-galactosamine. Vox Sang. 1987;52(1-2):125-8.

Laboratory Reagent For Research Use Only.

JK 07/06/16

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com