



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>  
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)



Delivering You Today's Innovations for  
the Science of Tomorrow™

**Purified Rabbit anti-Mycobacterium tuberculosis antigen EsaT-6 (Rv3875)  
Polyclonal Antibody**

**CLX305AP**

**Lot:**

**Clone:** Polyclonal

**Isotype:** Rabbit None

**Specificity:** The polyclonal antibody reacts with EsaT-6 (6 kDa early secretory antigen target ESXA), an 11.8 kDa antigen encoded by gene Rv3875 of *Mycobacterium tuberculosis*. The antibody recognizes EsaT-6 of *Mycobacterium tuberculosis H37Rv* and *Mycobacterium bovis*; it reacts also with recombinant antigen produced in *Escherichia coli*.

**Immunogen:** EsaT-6 (Rv3875) of *Mycobacterium tuberculosis*

**Application:** **Western Blotting**  
*Recommended dilution:* 0.2 - 1 µg/ml  
*Application note:* Reducing conditions.

**Purity:** > 95% (by SDS-PAGE)

**Purification:** Purified from rabbit serum by protein-A affinity chromatography.

**Concentration:** 1 mg/ml

**Storage Buffer:** Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4

**Storage/Stability:** Store at 2-8°C. Do not use after expiration date stamped on vial label. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.

**References:** Unpublished.

*This product is made available in collaboration with Exbio Praha, a.s.*

**Laboratory Reagent For Research Use Only**

JV 06/28/10

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup>



[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)