



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Purified Mouse Anti-HLA-DR Monoclonal Antibody

**CLX470AP**

<b>Size:</b>	0.1 mg
<b>Clone:</b>	L243
<b>Isotype:</b>	Mouse IgG <sub>2a</sub> kappa
<b>Specificity:</b>	The mouse monoclonal antibody L243 recognizes specifically an extracellular epitope on HLA-DR molecules, both peptide-loaded and empty.
<b>Immunogen:</b>	Human B lymphocytes
<b>Species Reactivity:</b>	Dog, Non-human primates, Human
<b>Application:</b>	FC (QC tested), IP, WB, IHC(P), IHC(F), ICC
<b>Application Details:</b>	Flow cytometry: Recommended dilution: 1-4 µg/ml.
<b>Purification:</b>	Purified from cell culture supernatant by protein-A affinity chromatography.
<b>Concentration:</b>	1 mg/ml
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4.
<b>Storage / Stability:</b>	Store at 2-8°C. Do not freeze.

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2015 and ISO 13485:2016  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)

**References:**

Coral S, Pucillo C, Leonardi A, Fonsatti E, Altomonte M, Maio M: Triggering of HLA-DR antigens differentially modulates tumor necrosis factor alpha release by B cells at distinct stage of maturation. *Cell Growth Differ.* 1997 May;8(5):581-8.

Cantin R, Martin G, Tremblay MJ: A novel virus capture assay reveals a differential acquisition of host HLA-DR by clinical isolates of human immunodeficiency virus type 1 expanded in primary human cells depending on the nature of producing cells and the donor source. *J Gen Virol.* 2001 Dec;82(Pt 12):2979-87.

De Gassart A, Camosseto V, Thibodeau J, Ceppi M, Catalan N, Pierre P, Gatti E: MHC class II stabilization at the surface of human dendritic cells is the result of maturation-dependent MARCH I down-regulation. *Proc Natl Acad Sci U S A.* 2008 Mar 4;105(9):3491-6.

Kalka-Moll WM, Tzianabos AO, Bryant PW, Niemeyer M, Ploegh HL, Kasper DL: Zwitterionic polysaccharides stimulate T cells by MHC class II-dependent interactions. *J Immunol.* 2002 Dec 1;169(11):6149-53.

Muczynski KA, Ekle DM, Coder DM, Anderson SK: Normal human kidney HLA-DR-expressing renal microvascular endothelial cells: characterization, isolation, and regulation of MHC class II expression. *J Am Soc Nephrol.* 2003 May;14(5):1336-48.

Swiatek-de Lange M, Rist W, Stahl HF, Weith A, Lenter MC: Comment on "MHC class II expression identifies functionally distinct human regulatory T cells". *J Immunol.* 2008 Mar 15;180(6):3625; author reply 3626.

Ivanov A, Beers SA, Walshe CA, Honeychurch J, Alduaij W, Cox KL, Potter KN, Murray S, Chan CH, Klymenko T, Erenpreisa J, Glennie MJ, Illidge TM, Cragg MS: Monoclonal antibodies directed to CD20 and HLA-DR can elicit homotypic adhesion followed by lysosome-mediated cell death in human lymphoma and leukemia cells. *J Clin Invest.* 2009 Aug;119(8):2143-59.

Bouillon M, El Fakhry Y, Girouard J, Khalil H, Thibodeau J, Mourad W: Lipid raft-dependent and -independent signaling through HLA-DR molecules. *J Biol Chem.* 2003 Feb 28;278(9):7099-107.

**Laboratory Reagent For Research Use Only**