



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Axygen® HybridRack™ Pipet Tips

Improve your environmental footprint with Corning's innovative hybrid design



Have you ever wanted your lab to use less plastic? Have you been looking for an easy way to be more sustainable around the lab?

Just like you, Corning takes sustainability seriously and is always looking for ways to reduce our carbon footprint. The HybridRack does just that, marrying functionality and sustainability by pairing a polypropylene deck with a paperboard base to allow for more sustainability in science. All without sacrificing quality, comfort, or reliability.

Innovative paperboard and plastic hybrid design uses


70% less plastic

than a traditional tip rack

Take care of the environment with Corning

Paperboard base can be easily recycled using most curbside recycling programs which can significantly reduce the amount of plastic being thrown away in the lab.

- ▶ Sterile
- ▶ Nonpyrogenic
- ▶ RNase-/DNase-free
- ▶ Human gDNA- and PCR-inhibition-free
- ▶ Available with polyethylene filter
- ▶ Available with Axygen Maxymum Recovery® surface to improve accuracy and ensure maximum sample recovery



Functionality meets sustainability

The unique hybrid design creates a stable platform allowing end users to pipet directly from the paperboard base – even with a multi-channel instrument.



This product is included in Corning EcoChoice™. All products included in the program are produced, packaged, or distributed in an environmentally friendly manner following United States FTC Guidelines. To learn more, go to www.corning.com/EcoChoice.

- ▶ Same consistent and high quality performance you have come to expect from Axygen pipet tips
- ▶ Unique design allows for use as a filter tip reload system
- ▶ Plastic deck is compatible with Axygen MultiRack and Corning® DeckWorks™ pipet tip racks
- ▶ Made in the USA in an ISO 13485 and ISO 14001 registered facility

Ordering Information

Products may not be available in all markets.

Axygen® HybridRack™ Pipet Tips

Racked, Filtered, Axygen Maxymum Recovery® surface

Cat. No.	Description	Quantity	Similar Plastic Rack Cat. No.
ERF-10XT-L-R-S	10 µL HybridRack extended length filtered tip, racked, Maxymum Recovery surface, sterile	960/pk, 5 pks/cs	MRF-10XT-L-R-S
ERF-20-L-R-S	20 µL HybridRack filtered tip, racked, Maxymum Recovery surface, sterile	960/pk, 5 pks/cs	MRF-20-L-R-S
ERF-100-L-R-S	100 µL HybridRack filtered tip, racked, Maxymum Recovery surface, sterile	960/pk, 5 pks/cs	MRF-100-L-R-S
ERF-200NX-L-R-S	200 µL HybridRack filtered tip, racked, Maxymum Recovery surface, sterile	960/pk, 5 pks/cs	MRF-200NX-L-R-S
ERF-300XT-L-R-S	300 µL HybridRack extended length filtered tip, racked, low binding, Maxymum Recovery surface, sterile	960/pk, 5 pks/cs	MRF-300XT-L-R-S
ERF-1000XT-L-R-S	1000 µL HybridRack extended length filtered tip, racked, Maxymum Recovery surface, sterile	960/pk, 5 pks/cs	MRF-1000XT-L-R-S

Racked, Filtered

Cat. No.	Description	Quantity	Similar Plastic Rack Cat. No.
ERF-10XT-R-S	10 µL HybridRack extended length filtered tip, racked, sterile	960/pk, 5 pks/cs	MRF-10XT-R-S
ERF-20-R-S	20 µL HybridRack filtered tip, racked, sterile	960/pk, 5 pks/cs	MRF-20-R-S
ERF-100-R-S	100 µL HybridRack filtered tip, racked, sterile	960/pk, 5 pks/cs	MRF-100-R-S
ERF-200NX-R-S	200 µL HybridRack filtered tip, racked, sterile	960/pk, 5 pks/cs	MRF-200NX-R-S
ERF-300XT-R-S	300 µL HybridRack extended length filtered tip, racked, sterile	960/pk, 5 pks/cs	MRF-300XT-R-S
ERF-1000XT-R-S	1000 µL HybridRack extended length filtered tip, racked, sterile	960/pk, 5 pks/cs	MRF-1000XT-R-S

Racked, Non-filtered

Cat. No.	Description	Quantity	Similar Plastic Rack Cat. No.
ER-10XT-R-S	10 µL HybridRack extended length tip, racked, sterile	960/pk, 5 pks/cs	MR-10XT-R-S
ER-200-R-S	200 µL HybridRack tip, racked, sterile	960/pk, 5 pks/cs	MR-200-R-S
ER-300NX-R-S	300 µL HybridRack tip, racked, sterile	960/pk, 5 pks/cs	MR-300-R-S
ER-300XT-R-S	300 µL HybridRack extended length tip, racked, sterile	960/pk, 5 pks/cs	MR-300XT-R-S
ER-1000XT-R-S	1000 µL HybridRack extended length tip, racked, sterile	960/pk, 5 pks/cs	MR-1000XT-R-S

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110.
Outside the United States, call +1.978.442.2200.