



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## NCI-H209 [H209] Cells Complete Medium

Cat. No: CM-0393  
Size: 125mL×4

### General Information

|                      |   |
|----------------------|---|
| Product Form         | Liquid  |
| Concentration        | Ready-to-use  |
| Component            | RPMI-1640( <a href="#">PM150110</a> )+10% FBS+1%<br>P/S( <a href="#">PB180120</a> ) |
| Bacterial detection  | Negative  |
| Fungal detection     | Negative  |
| Mycoplasma detection | Negative  |
| Endotoxin content    | ≤5  |
| Storage              | 2-8°C, Shading light  |
| Shipping             | Ice bag   |
| Expiration date      | 3 months  |

### Notes

1. This product is for research use only.
2. This product is sterilized by 0.22 μ m filtration.
3. It is necessary to pay attention to the aseptic operation and avoid the pollution during the culture.
4. It is not suitable for long time storage at room temperature or 4 °C.
5. This product contains serum, penicillin and streptomycin. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.