



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Penicillin-Streptomycin-Neomycin Solution (PSN) (100×)

Cat. No: PB180123
Size: 100mL

General Information

Product Form	Liquid
Concentration	100×
pH	6.1-6.4
Neomycin sulfate	10mg/ml
Dissolvent	10mM PBS (pH6.2)
Streptomycin sulfate	5mg/ml
Penicillin G	5kU/ml
Antimicrobial spectrum	Gram-positive bacteria & Gram-negative bacteria
Storage	-20°C, Shading light
Shipping	Ice bag or dry ice
Expiration date	12 months

Background

Penicillin-streptomycin solution mixture is the most commonly used antibiotic to prevent microbial contamination in vitro. Penicillin can interfere with the synthesis of bacterial cell wall, especially for Gram-positive bacteria. Streptomycin could bind to ribosomal ribosome for 30S and inhibit the synthesis of bacterial protein. It was effective for Gram-negative bacteria and Gram-positive bacteria, but especially for Gram-negative bacteria. Neomycin can bind to the 30S and 50S subunits of the bacterial ribosome to cause protein miscoding, and is effective against both Gram-negative and Gram-positive bacteria. Penicillin solution is sensitive to temperature and pH, easy to degrade at room temperature, and needs cryopreservation, the stability of penicillin solution is the most stable when pH is 6.0 - 6.5. Streptomycin was relatively stable, and pH 5.0-7.5 was the most stable, while neomycin sulfate is more stable.

Notes

1. This product is for research use only.
2. This product is sterilized by 0.22 μ m filtration.
3. It is necessary to pay attention to the aseptic operation and avoid the pollution.
4. Avoid repeated freezing and thawing. It is recommended to defrost at 2-8°C and aliquot the solution and store at -20°C.
5. This product is a concentrated solution. Please dilute it as needed.