



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

0.02% EDTA Solution (Versene Solution)

Cat. No: PB180320
Size: 100mL

General Information

Product Form	Liquid
Concentration	0.02%
pH	7.2-7.4
Phenol red	negative
Dissolvent	10mM PBS
EDTA	0.2g/L
Storage	2-8°C
Shipping	RT
Expiration date	18 months

Background

EDTA is a chemical chelating agent that can chelate divalent metal ions such as Ca^{2+} and Mg^{2+} . Its solution is also called Versene solution. Because Ca^{2+} and Mg^{2+} are important components of cell membranes and have the function of agglutinating cells, EDTA can destroy the connection between cells and promote the dissociation of cells into single cells. EDTA is less toxic, cheap, convenient to use, and does not damage cell surface molecules. It only binds to Ca^{2+} and Mg^{2+} , so it is suitable for cell digestion where cell surface molecules need to be detected. It is also often mixed with trypsin solution in a certain proportion.

Notes

1. This product is for research use only.
2. This product is sterilized by 0.22 μm filtration.
3. It is necessary to pay attention to the aseptic operation and avoid the pollution.
4. EDTA will not be neutralized, and the digested cells should be washed with PBS, otherwise the cells will be hard to adhere.