



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Krebs-Ringer Solution

Cat. No: PB180347
Size: 500mL

General Information

Product Form	Liquid
Concentration	1×
pH	7.2-7.4
D-glucose	1400 mg/L
NaHCO ₃	1370mg/L
Phenol red	Negative
CaCl ₂ (anhydrous)	370mg/L
MgCl ₂ •6H ₂ O	470mg/L
Storage	2-30°C
Shipping	RT
Expiration date	36 months

Background

Balanced Salt Solution (Physiological Solution) have the properties of buffer solution (regulate pH), normal saline (maintain osmotic pressure) and culture medium (provide nutrition). It can meet the basic needs of survival and metabolism of tissues, organs or cells in vitro. A small amount of phenolic red was added to some equilibrium salt solutions to indicate the pH change of the solution. Krebs-Ringer solution is a kind of balanced salt solution in biochemistry. The main components are NaCl, KCl, NaH₂PO₄, MgCl₂, CaCl₂, NaHCO₃ and glucose, is suitable for various mammalian tissues, such as nerve tissue, isolated heart perfusion, vascular or intestinal perfusion, etc.

Notes

1. This product is for research use only.
2. This product is sterilized by 0.22 μ m filtration.
3. It is necessary to pay attention to the aseptic operation and avoid the pollution.
4. Since the solution contains calcium and magnesium, there maybe precipitation produced, which will be dissolved after 37°C water bath incubation.