



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Leibovitz's L-15 Complete Medium

Cat. No: PM151010C  
Size: 125mL×4

### General Information

Product Form	Liquid
Concentration	Ready-to-use
Component	Leibovitz's L-15( <a href="#">PM151010</a> )+20% FBS+1% PS( <a href="#">PB180120</a> )+2mM L-Alanyl-L-Glutamine( <a href="#">PB180419</a> )
Bacterial detection	Negative
Fungal detection	Negative
Mycoplasma detection	Negative
Endotoxin content	≤5
Storage	2-8°C, Shading light
Shipping	Ice bag
Expiration date	3 months

### Notes

1. This product is for research use only.
2. This product is sterilized by 0.22 μ m filtration.
3. It is necessary to pay attention to the aseptic operation and avoid the pollution during the culture.
4. It is not suitable for long time storage at room temperature or 4 °C.
5. This product contains serum, penicillin and streptomycin. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.