



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Exalpha Biologicals, Inc.

Human Heart Tissue Control Lysate

CAT. No.
X1634C
SIZE

100 μ l

CONCENTRATION

1 mg/ml

APPLICATIONS

Western Blot

COMMENTS

Human heart tissue lysate is ready to load on SDS-PAGE for Western blotting.

FORMULATION

Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol

STORAGE CUSTOMER

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

STABILITY

Products are stable for one year from purchase when stored properly

USE

Use 20 μ l of human heart tissue lysate per lane. Lysate should be boiled prior to loading to ensure total denaturation of lysate. Boil for 10 minutes prior to loading.

For research use only. Not for use in human diagnostics or therapeutics.

Exalpha Biologicals, Inc. 5 Clock Tower Place, Suite 255 Maynard, MA 01754
Tel: 800.395.1137 Fax: 866.924.5100 www.exalpha.com info@exalpha.com

LAST MODIFIED
4/4/2007