



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Exalpha Biologicals, Inc.

## Human Placenta Tissue Control Lysate

**CAT. No.**  
X1635C  
**SIZE**

100  $\mu$ l

**CONCENTRATION**

1 mg/ml

**APPLICATIONS**

Western Blot

**COMMENTS**

Human placenta tissue lysate is ready to load on SDS-PAGE for Western blotting.

**FORMULATION**

Provided in 10% glycerol, 0.063 M Tris-HCl (pH 6.8), 2% SDS and 0.002% bromophenol blue, 5% 2-mercaptoethanol

**STORAGE CUSTOMER**

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

**STABILITY**

Products are stable for one year from purchase when stored properly

**USE**

Use 20  $\mu$ l of human placenta tissue lysate per lane. Lysate should be boiled prior to loading to ensure total denaturation of lysate. Boil for 10 minutes prior to loading.

**For research use only. Not for use in human diagnostics or therapeutics.**

Exalpha Biologicals, Inc. 5 Clock Tower Place, Suite 255 Maynard, MA 01754  
Tel: 800.395.1137 Fax: 866.924.5100 www.exalpha.com info@exalpha.com

**LAST MODIFIED**  
4/4/2007