



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Chicken anti Bovine Serum Albumin (BSA), conjugated with Agarose

[exalpha.com/products/chicken-anti-bovine-serum-albumin-bsa-conjugated-with-agarose/ABSA-AGA-1](https://www.exalpha.com/products/chicken-anti-bovine-serum-albumin-bsa-conjugated-with-agarose/ABSA-AGA-1)

Catalog number: **ABSA-AGA-1**

Product Type	Primary antibodies
Units	1.0 ml bed volume
Host	Chicken

Source

Chickens were immunized with highly purified bovine serum albumin isolated from normal pooled bovine serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbant and conjugated to N-hydroxysuccinimide (NHS)-activated agarose beads. Beads were washed extensively after the blocking of the residual NHS sites.

Product

Gel; 1.0 ml bed volume

Product Form: Immunoprecipitation Gel

Formulation: Phosphate buffered saline, pH 7.2 with 0.075% sodium azide

Purification Method: Egg yolk-derived purified IgY (non affinity purified)

Concentration: N/A

Applications

This product is suitable for immunoprecipitation and immunodepletion procedures. Optimal working dilutions should be determined for your particular assay conditions.

Storage

Product should be stored at 4-8°C. DO NOT FREEZE

Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional

information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpha Biologicals accepts no liability for any inaccuracies or omissions in this information.

Safety Datasheet(s) for this product:

EA_Sodium Azide