



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Chicken anti Botulinum Neurotoxin Type E Light Chain, conjugated with Biotin

[exalpa.com/products/chicken-anti-botulinum-neurotoxin-type-e-light-chain-conjugated-with-biotin/ACbNELC-B](https://www.exalpa.com/products/chicken-anti-botulinum-neurotoxin-type-e-light-chain-conjugated-with-biotin/ACbNELC-B)

Catalog number: **ACbNELC-B**

Product Type	Primary antibodies
Units	0.1 mg
Host	Chicken
Application	ELISA Western Blotting

Source

Chickens were immunized with recombinant Botulinum Neurotoxin Type E Light Chain (Clostridium botulinum). After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY purified. Anti-Botulinum Neurotoxin Type E Light Chain was immunoaffinity purified using the immunogen immobilized on a solid support and biotinylated.

Product

Affinity purified antiserum in phosphate buffered saline, pH 7.2, 50% glycerol with 0.1% 2-chloroacetamide.

Product Form: Biotin conjugated

Purification Method: Egg-yolk derived IgY Affinity Purified

Concentration: See vial for concentration.

Applications

Optimal working dilutions should be determined for your particular assay conditions.

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles.

Shipping Conditions: Ship at ambient temperature.

Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may

contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpha Biologicals accepts no liability for any inaccuracies or omissions in this information.

Safety Datasheet(s) for this product:

EA_Thimerosal EA_Glycerol