



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Chicken anti Cervid IgG (H+L), conjugated with FITC

[exalpa.com/products/chicken-anti-cervid-igg-hl-conjugated-with-fitc/ACerIG-F](https://www.exalpa.com/products/chicken-anti-cervid-igg-hl-conjugated-with-fitc/ACerIG-F)

Catalog number: **ACerIG-F**

Product Type	Secondary antibodies
Units	0.5 mg
Host	Chicken
Species Reactivity	Cervid
Application	Immunohistochemistry

Source

Chickens were immunized with highly purified cervid IgG isolated from deer and elk serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbant and conjugated to fluorescein isothiocyanate.

Product

Product Form: Fluorescein conjugated

Formulation: Phosphate Buffered Saline, 50% Glycerol with 0.075% sodium azide

Purification Method: Egg-yolk derived IgY Affinity Purified

Concentration: See vial for concentration.

Applications

Optimal working dilutions should be determined for your particular assay conditions.

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles.

Shipping Conditions: Ship at ambient temperature

Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional

information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpa Biologicals accepts no liability for any inaccuracies or omissions in this information.

References

1. Duscher, GG, Wetscher, M, Baumgartner, R and Walder G. 2015. Roe deer sera used for TBE surveillance in Austria. Ticks Tick Borne Dis. 6(4):489-93.

Safety Datasheet(s) for this product:

EA_Sodium Azide EA_Glycerol