



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



## Chicken anti Human IgG (H/L), conjugated with Biotin

[exalpha.com/products/chicken-anti-human-igg-h-l-conjugated-with-biotin/AHIgG-B](https://www.exalpha.com/products/chicken-anti-human-igg-h-l-conjugated-with-biotin/AHIgG-B)

Catalog number: **AHIgG-B**

Product Type	Secondary antibodies
--------------	----------------------

Units	0.5 mg
-------	--------

Host	Chicken
------	---------

Application	ELISA Western Blotting
-------------	---------------------------

### Source

Chickens were immunized with highly purified normal human IgG. After multiple immunizations, eggs were collected and the IgY fraction prepared. This fraction was, immunoaffinity purified from an Immunoabsorbant and conjugated to biotin.

### Product

*Product Form:* Biotin conjugated

*Formulation:* Phosphate buffered saline, pH 7.2, 50% glycerol with 0.01% Thimerosal

*Purification Method:* Egg-yolk derived IgY Affinity Purified

*Concentration:* See vial for concentration.

### Applications

Optimal working dilutions should be determined for your particular assay conditions.

### Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles.

### Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpha Biologicals accepts no liability for any inaccuracies or omissions in this information.

**References**

1. Gorkin III, R., Nwankire, C.E., Gaughran, J., Zhang, X. Donohoe, G.G., Rook, M., Kennedy R. and J. Ducreé. 2012. Centrifugo-pneumatic valving utilizing dissolvable films. *Lab Chip* 12:2894-2902.
2. Nwankire, C.E., Donohoe, G.G., Zhang, X., Siegrist, J., Somers, M., Kurzbuch, D., Monaghan, R., Kitsara, M., Burger, R., Hearty, S., Murrell, J., Martin, C., Rook, M., Barrett, L., Daniels S., McDonagh, C., O'Kennedy, R., and J. Ducree. 2013. At-line bioprocess monitoring by immunoassay with rotationally controlled serial siphoning and integrated supercritical angle fluorescence optics. *Anal. Chim. Acta* 781:54-62.

**Safety Datasheet(s) for this product:**

EA\_Thimerosal EA\_Glycerol