



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Chicken anti Sheep IgG (H+L)

[exalpha.com/products/chicken-anti-sheep-igg-hl/ASIG](https://www.exalpha.com/products/chicken-anti-sheep-igg-hl/ASIG)

Catalog number: **ASIG**

Product Type	Secondary antibodies
--------------	----------------------

Units	0.5 mg
-------	--------

Host	Chicken
------	---------

Application	ELISA Western Blotting
-------------	---------------------------

Source

Chickens were immunized with highly purified sheep IgG isolated from normal pooled sheep serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbent.

Product

Product Form: Affinity purified

Formulation: Phosphate buffered saline, pH 7.2 with 0.075% sodium azide

Purification Method: Egg-yolk derived IgY Affinity Purified

Concentration: See vial for concentration.

Applications

This product is suitable for ELISA and blotting procedures. Optimal working dilutions should be determined for your particular assay conditions. Storage

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles.

Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional

information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Exalpa Biologicals accepts no liability for any inaccuracies or omissions in this information.

References

1. Rubenstein, R., Chang, B., Gray, P., Piltch, M., Bulgin, MS., Sorenson-Melson S. and MW Miller. 2011. Prion Disease Detection, PMCA Kinetics, and IgG in Urine from Sheep Naturally/Experimentally Infected with Scrapie and Deer with Preclinical/Clinical Chronic Wasting Disease. *J. Virol.* 85(17):9031-8. Epub 2011 Jun 29.

Safety Datasheet(s) for this product:

EA_Sodium Azide