



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

10X Citrate buffer pH8.0

Cat No. GTX03595

Application IHC-P

Package
500 ml

PRODUCT

Summary

Shi et al. (1991) have reported that treating tissue sections with a heavy metal solution in a microwave oven can cover masked antigens significantly. However, heavy metals in the solution increase the risk of toxic exposure to lab personnel. To reduce the risk, we have developed an antigen unmasking solution which is free of heavy metals. The reagent is provided in a convenient 10X concentrate. Use of this antigen recovery buffer avoids unnecessary heavy metal exposure to lab personal and handling and disposal issues.

PROPERTIES

Form Liquid

Storage Store at room temperature.

Conjugation Unconjugated

Note For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).