



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Multiple Tags Positive Control (6xHis-MBP-DDDDK-V5-HSV-HA-T7-Myc-S)


Cat. No. GTX130342-pro

Package
100 µg, 25 µg

PROPERTIES

Form	Liquid
Buffer	1X SDS-PAGE sample loading buffer
Preservative	No Preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Expression System	E. coli
Purity	> 90%
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX130342-pro Image

Lane1: 1 µg of negative control
 Lane2: 1 µg of multiple tags positive control (GTX130342-pro) blots were probed with various GeneTex's Tag antibodies and then probed with goat polyclonal secondary antibody to rabbit IgG - H&L (HRP) (GTX213110-01) or goat polyclonal secondary antibody to mouse IgG - H&L (HRP) (GTX213111-01).



For full product information, images and publications, please visit our [website](#).