



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# AP Enhancer

**Cat No. GTX20671**

**Application** IHC-P, IHC-Fr

**Package**  
250 ml

## PRODUCT

**Product Description** Alkaline Phosphatase Enhancer

## APPLICATION

### Application Note

i) Once tissue sections have been incubated with streptavidin-alkaline phosphatase, wash them with buffer thoroughly. ii) Wipe the glass to remove excess of buffer and add enough drops of the alkaline phosphatase enhancer solution to cover the tissue sections. Alternately, pour alkaline phosphatase enhancer solution in a coplin jar and dip the slides for 1 minute. iii) Incubate for 1 minute at room temperature. iv) Drain excess alkaline phosphatase enhancer and add substrate/chromogen onto the tissue sections without any wash between the two steps. Precautions: Avoid contact with clothes and exposed skin. If accidentally contacted, flush with tap water immediately. Do not take internally.

## PROPERTIES

**Form** Liquid

**Buffer** 250 ml of stable, ready to use solution.

**Storage** Store at room temperature.

**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).