



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# DAB Enhancer

**Cat No. GTX20675**

**Application** IHC-P, IHC-Fr

**Package**  
10 ml

## PRODUCT

**Product Description** DAB Enhancer

## APPLICATION

### Application Note

i) Add one drop (50µl) of DAB enhancer in 1 ml of DAB incubation buffer. ii) Add DAB chromogen and hydrogen peroxide to the above buffer. iii) Incubate tissue sections with the incubation medium.

## PROPERTIES

**Form** Liquid

**Buffer** 10 ml concentrated solution.

**Storage** Store at room temperature.

**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).