



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TissueBond slide adhesive

Cat No. GTX29909

Application	IHC-P, IHC-Fr
--------------------	---------------

Package	7 ml
----------------	------

PRODUCT

Product Description	TissueBond slide adhesive
----------------------------	---------------------------

APPLICATION

Application Note

Prepare of working solution:

- i) Working solution should be prepared just prior to use.
- ii) Add the entire content of TissueBond in a Acetone compatible container.
- iii) Add dry acetone to make a volume of 350 ml and stir well. This solution is enough for coating 500 slides.

*Once the bottle has been opened, use it immediately. Prolonged storage is not recommended. Discard unused working solution.

Application Note:

- 1) Place slides in a metallic slide holder. Clean slides by immersing them in acetone for 1 minute.
- 2) Immerse slides in TissueBond working solution for 30 sec.
- 3) Remove and drain the slides. Rinse the slides thoroughly with deionized water three times. Change the water after every rinse. Avoid creating water bubbles by gentle agitation of slide holder before drying.
- 4) Air dry slides at room temperature. Once dried, slides can be used immediately or can be stored in a box at room temperature for future use.

PROPERTIES

Form	Liquid
-------------	--------

Storage	Keep as concentrated solution, aliquot and store at 4°C.
----------------	--

Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	--



For full product information, images and publications, please visit our [website](#).