



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# ImmunoHistoMount™

**Cat No. GTX30922**

<b>Application</b>	IHC
--------------------	-----

**Package**  
30 ml

## PRODUCT

<b>Product Description</b>	ImmunoHistoMount™
----------------------------	-------------------

## APPLICATION

### Application Note

Apply 3-4 drops of mounting medium directly on top of the specimen and spread out evenly by tilting the slide back and forth or spread evenly with a 0.2 ml plastic pipette tip making sure the tissue is not touched. Excess medium can be removed by touching the edges of slide against a paper towel. Let stand at room temperature for about 1 to 2 hours. Or the slides can be heated at 37-40oC for 40-60 minutes. The slides are ready for visualization under a microscope.

## PROPERTIES

<b>Form</b>	Liquid
-------------	--------

<b>Storage</b>	This product is stored at room temp, however for long term storage 2-8oC is recommended. The pH of this product is more stable at 2-8oC. DO NOT FREEZE.
----------------	---

<b>Note</b>	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	--



For full product information, images and publications, please visit our [website](#).