



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Organo/Limonene Mount™

Cat No. GTX30923

Application IHC

Package
30 ml

PRODUCT

Product Description Organo/Limonene Mount™

APPLICATION

Application Note

Dehydrate tissue or cell smear slides with dehydrating agents e.g. 30 to 100% alcohol followed by xylene, toluene or any other dehydrating reagents. Add 2-3 drops of Organo mounting medium. Apply coverslip carefully, avoiding air bubbles. Visualize slide under a microscope. The slides can be dried by leaving at room temperature overnight or by heating at 37-40oC for several hours.

PROPERTIES

Form Liquid

Storage Room temperature

Note For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).