



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Alkaline Phosphatase Enhancer

**Cat No. GTX30944**

**Application** ICC/IF, IHC, ISH

**Package**  
250 ml

## PRODUCT

**Product Description** Alkaline Phosphatase Enhancer

## APPLICATION

### Application Note

1. After the alkaline phosphatase reaction, wash slide with buffer, and distilled/ deionized water.
2. Add few drops of this ready to use Enhancing Reagent to cover the tissue/cell smear, for 1-2 minutes and discard. DO NOT WASH SLIDE.
3. Add few drops of ready to use alkaline phosphatase chromogen, follow instructions as outlined for these chromogens.

## PROPERTIES

**Storage** Room temperature

**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).