



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Trident Universal Protein Blocking Reagent (animal serum free)

**Cat No. GTX30963**

**Application** WB, ICC/IF, IHC-P, ELISA, IHC

**Reference ( 2 )**

**Package**

100 ml

## PRODUCT

**Product Description** Protein Blocking Reagent (animal serum free)

## APPLICATION

### Application Note

Ready to use. Western blot: Immerse the filter in a sufficient amount of Blocking Reagent and incubate for 30 minutes at RT with gentle shaking. IHC/ICC: Incubate tissue sections/ cell smears in Blocking Reagent for 10 minutes at RT or 37°C. Note: Optimal incubation times should be determined by the researcher.

## PROPERTIES

**Form** Liquid

**Buffer** 20mM Tris 150mM, sodium chloride, 0.02% Triton X100, pH 7.4, and 0.5% Proclin 300 was added as a preservative. It does not contain SDS or sodium azide or any other preservative.

**Storage** 2-8°C

**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).